

MRI Compatibility Testing for Accessories

Sound Imaging Inc's MRI Comfort Suite Items were tested at Sharps & Children's MRI on a 750 3T GE MRI and on a 750W 3T GE MRI at Scripps Memorial Hospital for Compatibility in a 3T MRI System. Daniel Rubalcaba, with the guidance of Bob Bell Ph.D., was present for the testing.

Items were checked to see if they could be used near the bore of the MRI Magnet. The units were measured to see if there was any pull, and tested for acceptable distance and use in the 3T MR Environment.

Test dates: 04/30/19 & 05/01/19

The following list describes the Items and the results:

Model Number	Description	Result	Part of:
SI-MRISTE-TND	MRI Stereo Transducer	MRI Conditional for 3T	MRI Stereo System
SI-SAMMCAM-1	MRI Camera w/ Lens	MRI Conditional for 3T	SAMM System
SI-SAMMCAW-1	MRI Camera White	MRI Conditional for 3T	SAMM White
SI-SAMMCAM-2	MRI Camera w/ Lens	MRI Conditional for 3T	SAMMmd

Chart Revised 05/02/2019

SI-03-26-00 iso doc

MRITesting.doc

Regards,



Daniel Rubalcaba

Director of Technology